crystalign

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Splints made from optimill crystalign

optimill crystalign is certified for use in the patient's mouth for orthodontic splints. The milling blanks have a very high material homogeneity due to the industrial polymerization process. The combination of breaking strength and elasticity, which is well-balanced for clinical use, guarantees excellent long-term stability. optimill crystalign has admirable abrasion resistance. The material ensures high process reliability, no mixing errors, no polymerization shrinkage and no unpleasant smell.

Indications for optimill crystalign

optimill crystalign is intended for use in the patient's mouth for the creation of orthodontic splints, which are designed by a dental technician and manufactured using CAD-CAM technology.

Contraindications for optimill crystalign

Do not use this product in sensitive or allergic patients. NOTE: Please pass all designated information to the dentist, so that the patients can be consulted according to possible sensitivities or allergies prior to the treatment. The material should not be used for any other purposes than the generative CAD / CAM production of orthodontic splints.

Side effects for optimill crystalign

Correct processing and application by the dentist and dental technician observed no side effects. Immune reactions (for example allergies) cannot be completely ruled out. Please inform us of about occuring side effects.

Processing

To process optimill crystalign use only carbide burs. Please note that the improper use of milling tools, polishing brushes, steam jets and water baths can cause the material to overheat. This may result in damage of the material.

Milling

Use fine-toothed tungsten carbide bur to separate the milled restoration from the disc. Cut the retaining bars between the product and disc carefully and without applying pressure. Orthodontic splints can be pre polished with suitable silicone polishes and goat hair brushes. To achieve a high shine polish, it is recommended, to use the diamond polishing pastes. In order to prevent plaque formation, a careful polishing is an essential prerequisite.

Cleaning

After preparation, clean the finished objects with compressed air or in an ultrasonic bath. Steam cleaners are not recommended as deformations of the product may occur due to the high temperature and pressure.

Properties

Chemical characterization	-	PMMA Polymethylmethacrylat / Pigment
Flexibility	MPa	≈ 78
Water absorption	μg/mm3	≈ 19-25
E modules	MPA	≈ 2988
Residual monomer content	%	0,02

Disposal

Orthodontic splints made with optimill crystalign and leftovers are insoluble in water, inactive, pose no threat to ground water and can thus be disposed of in normal domestic waste.

Symbols

See instruction for use

 $\overrightarrow{\uparrow}$

Keep in dry place

Lot number

Keep away from sunlight

REF Product number

Distributed by: dentona AG

Otto-Hahn-Straße 27 | 44227 Dortmund Phone: +49 (0) 231 - 55 56 - 0 Fax: +49 (0) 231 - 55 56 - 30

manufacturer: on dent TIBBI MALZ. DIS PRO. ITH. IHR. SAN, ve TIC LTD. STI. 35070 Bornova – IZMIR Turkey

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Vertrieb durch: dentona AG
Hersteller: on dent TIBBI MALZ. DIS PRO. ITH. IHR. SAN, ve TIC LTD. STI. 35070 Bornova – IZMIR Turkey