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### Instructions for use

#### dentona-tray LC

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Light curing resin in form of contoured wafers for the fabrication of bite  
plates and base plates for needle-point tracing

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
## Application

- The situation model should be prepared as required and insulated with an alginate insulation. Undercutting areas of the model should be blocked out before.
- Open the can at the handle tab, remove one plate incl. protective fleece and immediately reseal the can tightly.
- The protective fleece is removed from the sheet.
- The sheet is adapted to the insulated model and trimmed at the edges.
- The handle can be formed from the excess material. The handle is not used when making a base plate.
- The sheet is then polymerized on the model in a blue light or UVA light device (e. g. Lumos by dentona®) for 5 minutes. The instructions for use of the light unit used must be observed!
- The sheet is removed from the model and polymerized from the back for another 5 minutes.
- The individual impression tray made in this way can be processed immediately.
- The inhibition layer (adhesive layer on the top of the material) allows subsequent correction (e. g. reinforcement of the sheet). The inhibition layer can be easily removed with a moistened tissue with an alcoholic cleaning agent.

## Delivery units

| product         | version     | quantity          | ordernumber |
|-----------------|-------------|-------------------|-------------|
| dentona-tray LC | transparent | 1 can á 50 pieces | 53061       |
|                 | blue        | 1 can á 50 pieces | 53062       |
|                 | pink        | 1 can á 50 pieces | 53060       |
| dentona-base LC | pink        | 1 can á 50 pieces | 53070       |

## User profile

These are products for professional use by dental technicians or dentist. It is intended for use in the dental laboratory or dental practice. No special training is required. The finished products made of the material (individual impression trays, base plates, support pin registrations) are intended for professional use by dentists in a dental practice. Special training is not required.

## Patient group

The use of these products is not limited to any particular patient group. Contraindications must be observed.

## Undesirable side effects

The products contain methacrylates. These may cause allergic skin reactions. In rare cases sensitization is possible.

## Interactions

No interactions are known.

## Special storage instructions

Keep the can tightly closed. Store protected from light and at 4-25 °C.

## Shelf life

3 years from manufacture. Do not use after the expiry date stated on the label!

## Waste treatment

Do not allow product or primary packaging which has not been emptied of all residues to reach sewage system, water, subsoil or soil. Disposal according to official regulations. Non-contaminated and completely empty packaging as well as fully cured products can be disposed of with normal municipal waste (household waste). Observe waste separation.

## General safety instructions

 *Observe safety data sheet!*

H317: May cause an allergic skin reaction. Avoid contact of uncured material with skin/mucous membrane and eyes. It may be slightly irritating and in rare cases may cause sensitization to methacrylates.

Use protective equipment: Safety glasses, suitable protective gloves with protection function to Meth-acrylates (e.g. Polyethylene-gloves, 4h-gloves). Do not wear contaminated work clothing outside the workplace.

After skin contact: Wash immediately with plenty of soap and water. Medical treatment necessary if symptoms occur (irritation of the skin).

After eye contact: Immediately flush with running water for 10 to 15 minutes with the eyelid open and consult an ophthalmologist. When grinding the polymerized sheets, use a suitable dust extractor and follow the instructions for use of the abrasive manufacturer, especially with regard to dust protection. If necessary, use a self-contained breathing apparatus.

All serious incidents related to the product must be reported to the manufacturer and the competent authority of the member state where the user and/or patient is established. Please direct any feedback regarding the product to: [reclamation@megadenta.de](mailto:reclamation@megadenta.de)

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## dentona-tray LC and dentona-base LC

**General description**

**dentona-tray LC** and **dentona-base LC** are pre-shaped, modelable, light-curingcomposites for the production of individual impression trays in the dental laboratory or dental practice. They features low shrinkage, low water absorption and high dimensional stability.

Available variants

**dentona-tray LC:** plate thickness approx. 2.2 mm in the colors transparent, blue and pink

**dentona-base LC:** plate thickness approx. 1.4 mm in pink

**Composition**

Mixture of different acrylates / methacrylates, fumed silica, glass beads, photochemical initiators, vaseline, pigment

### Intended use

**Warnings**

The products are light-sensitive, so that polymerization can already be triggered by exposure to daylight. Therefore, the material should not be exposed to light unnecessarily (removal from the can only immediately before adaptation to the model, immediate sealing of the can after removal).

**Intended use**

**dentona-tray LC** is intended for dental use only:

- for the production of individual impression trays for precision impressions for accurately fitting restorations with the highest marginal accuracy, e. g.
  - functional impressions
  - situational impressions for unfavorable dental arch and jaw shapes and/or deep-set bands

- dentona-base LC** plates are for dental use only
  - for production of a occlusal base plate
  - for intraoral support pin registration

**Clinical benefit**

The clinical benefit of light-curing impression plates lies indirectly in the replacement of the dental apparatus and thus in a restoration of masticatory function, since precision and/or functional impressions are essential for the accurately fitting production of dentures.

**Indications**

**dentona-tray LC**

- Precision impressions for accurately fitting restorations with the highest marginal accuracy, such as:
  - functional impressions
  - Situation impressions for unfavorable dental arch and jaw shapes and/or deep-set bands

**dentona-base LC**

- Bite registrations
- Support pin registrations

**Contraindication**

In case of known hypersensitivity to methacrylates, skin contact should be avoided.

#### Deutsch

#### Deutsch

#### Deutsch

#### Englisch