For perfect fit and wear comfort

flexisplint®

dentona flexisplint®

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Chargennummer

Indication:

Dental resin for the production of all kinds of flexible dental therapeutic bite splints.

Contraindication:

The material should not be used for any other purpose than the production of dental bite splints. Dentona flexisplint® is not suitable for the manufacture of denture bases. Unwanted biological reactions (such as allergies against material components) can occur in very rare cases. At known incompatibilities or known allergies to the material components bite splints made from dentona flexisplint® must not be used.

Hazard warnings:

- The use of dentona flexisplint® can produce harmful dusts that can irritate the eyes, skin and airways. Always ensure that the extraction system at your workplace is working properly. (See also **safety data sheet** at www.dentona.de)
- The mixing ratio (**1 part powder / 0.75 parts liquid**) and the temperature during processing (**5-9 °C**) must be strictly observed, otherwise the resin cannot cure properly.
- Both components should be refrigerated for at least 15 minutes before mixing.
- To avoid uncontrollable chemical reactions, only use dosing cups and a spatula that don't come into contact with any other resin.

The optimized ratio for a dentona flexisplint® splint: 7,5 g liquid / 10 g powder

Pack pressing technique

1. Prepare the model in the normal way, measure and mark the tooth equator. As dentona flexisplint® is a flexible material, the splint can be arranged slightly above the tooth equator.
2. The splint to be created is waxed-up in the articulator taking into account the chewing movement (for example with pink wax sheet), and then bed the model into a cuvette. Now cool the required amount of powder in a refrigerator, (both components must be cooled before mixing).
3. Isolate all exposed plaster parts of the embedding with alginate isolation. As soon as the isolation is dry the counterpart can be produced from the modelled wax splint. Once the counterpart gypsum has completely set, the entire flask should be warmed up in a warm water bath. This makes it easier to open the two halves of the cuvette.
4. Remove the wax from the modelled wax splint, scald the model and the counterpart, water and isolate with alginate isolation.

Important: Adhere exactly

5. Remove the liquid from the refrigerator immediately prior to processing.
6. Mixing ratio: 1 part powder / 0.75 parts liquid
7. Measure the dentona flexisplint® material in separate dosing cups. Bring together in a rubber mixing cup, premix rapidly and bubble-free and then immediately fill into one side of the counterpart until the resin gets visible on the other side.
8. As soon as dentona flexisplint® has a pasty consistency after approx. 5 minutes, the resin is placed in the cuvette on the row of teeth and the cuvette top is attached. Any visible white dots in the mixed paste are production-related and will no longer be visible after polymerization.
9. Now quickly screw the cuvette under the press and then polymerize for at least 45 minutes in a water bath at 2 bar and 50 °C (maximum 53 °C!).
10. The cuvette must be completely cooled down before opening! This is ensured by placing the cuvette in a cold bath of water for approximately 15 minutes.
11. Now carefully lift the splint from the model and perform the usual finishing and polishing.

Casting technique

Prepare the model in the normal way, measure and mark the tooth equator. As dentona flexisplint® is a flexible material, the splint can be arranged slightly above the tooth equator.

1. The splint to be created is waxed-up in the articulator taking into account the chewing movement (for example with pink wax sheet), and then bed the model into a cuvette. Now cool the required amount of powder in a refrigerator, (both components must be cooled before mixing).
2. Now prepare a silicone counterpart (for example dentona® 1:1 gum) and allow to cure. Once the silicone counterpart is fully cured it can be carefully lifted off. Cut out an excess channel with a sharp scalpel on the tuber or both sides on the retromolar triangle.
3. Remove the wax from the modelled wax splint, scald the model and the counterpart, water and isolate the still warm plaster model with alginate isolation (the silicone counterpart is then isolated with a thin layer of Vaseline). Once the isolation is dry and the model is completely cooled down, (simply place it in the refrigerator), place the counterpart back on the model and secure.
4. Remove the liquid from the refrigerator immediately prior to processing.
5. Mixing ratio: 1 part powder / 0.75 parts liquid

Important: Adhere exactly

6. Measure the dentona flexisplint® material in separate dosing cups. Bring together in a rubber mixing cup, premix rapidly and bubble-free and then immediately fill into one side of the counterpart until the resin gets visible on the other side

7. Now polymerise the model with the filled silicone counterpart as quickly as possible - for 45 minutes - in a pressure pot at 2 bar and a minimum of 50 °C (maximum 53 °C!).
8. Then the workpiece is removed from the pressure pot and cooled in cold water before the finished splint is taken off.
9. Now carefully lift the splint from the model and perform the usual finishing and polishing.

Application tips:

1. Before removing from the model and before setting in the patient's mouth the splint must be made slightly flexible by exposing to lukewarm water.
2. For pre-polishing, use a leather buff with pumice at 3000 rpm. For the high-gloss polish, the cotton swab with a conventional high-gloss polish at 1500 rpm. is suitable.
3. Severe undercut areas and interdental spaces should be sufficiently blocked out so that when the splint is completed it can be lifted without damage from the model.
4. After completion, the minimum wall thickness of the splints must be 1,0 mm.

Composition:

Mono- and polymethacrylate, THFMA, stabilizer, radical initiator, catalyst

Storage temperature before opening: 5-30 °C

Notes:

Restorations are custom-made products in accordance with Directive 93/42/EEC. Store in clear water for at least 12 hours before delivery. Specify LOT No. for any operation that requires identification of the material.

Care / cleaning of the splints:

The cleaning of the objects before further processing should not be done by excessive evaporation due to the high heat development. An ultrasonic bath with a maximum temperature of approx. 35 °C and a soft toothbrush are more suitable. Cleaning solutions and denture cleaners are not recommended. Repairs can be done in the conventional manner. Before that, the area to be repaired needs to be roughened, moistened with a bit of "dentona flexisplint®" liquid and then processed with fresh dentona flexisplint® material.

Storage recommendation for the splint:

Avoid direct sunlight. When not in use, the finished splint should always be stored in clear water.

Disposal:

Dentona flexisplint® restorations are water-insoluble, inactive and do not contain any danger to the groundwater and can therefore be disposed of with normal household waste.

Use:

Use by qualified personnel only. Rx only.

Patient target group: adults and adolescents

Symbole

See instruction for use

Attention: harmful to health

Use by / Expiry date

Attention: health risk

CE mark with number of notified body

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