

optiprint® Eignungstabelle

Compatibility

Germany's favorite 3D Resin

		optiprint model	optiprint model align	optiprint gingiva 385	optiprint I-B-T 385	optiprint clara	optiprint guide	optiprint lumina	optiprint laviva	optiprint tray	optiprint zero	optiprint mock-up	optiprint match
<ul style="list-style-type: none"> ● Materialparameter vorhanden <i>Available material parameters</i> □ Materialparameter in Arbeit <i>Material parameters in progress</i> 													
Asiga Max (385 nm)		●	●	●	●	●	●	●	●	●	●	●	●
Asiga Pro 4K (385 nm)		●	●	●	●	●	●	●	●	●	●	●	●
EnvisionTEC D4K Pro (385 nm)		●	●	□	□	□	●	□	□	●	□	□	□
EnvisionTEC envisionOne (385 nm)		●	●	□	□	□	●	□	□	●	□	□	□
Ackuretta DENTIQ (405 nm)		●	●	●	●	●	●	●	●	●	●	●	●
Ackuretta SOL (405 nm)		●	●	●	●	●	●	●	●	●	●	●	●
Ackuretta FreeShape (405 nm)		●	●	●	●	●	●	●	●	●	●	●	●
HeyGears A2D (385 nm)		□	□	□	□	□	□	□	□	□	□	□	□
Microlay Eve (405 nm)		●	●	●	●	●	●	●	●	●	●	●	●
Microlay Versus (385 nm)		●	●	●	●	●	●	●	●	●	●	●	●
W2P SolLab 335/370/450 (385 nm)													
MiiCraft 125 (405 nm)		●		●	●					●			
MiiCraft Ultra (405 nm)													
Flashforge Hunter (405 nm)													
Phrozen Sonic XL 4K (405 nm)		●	●	●	●	●	●	●	●	●	●	□	□
Shining 3D Accu-Fab L4D (405 nm)		□	□	□	□	●	□	●	●	□	●	□	□
Shining 3D Accu-Fab D1s (405 nm)		□	□	□	□	●	□	●	●	□	●	□	□

dentona (04/2022)

Finden Sie regelmäßige updates unter www.dentona.de
 Visit www.dentona.de for regular updates