

optiprint® Eignungstabelle

Germany's favorite 3D Resin

	optiprint model	optiprint model align	optiprint gingiva 385	optiprint I-B-T 385	optiprint clara	optiprint guide	optiprint lumina	optiprint laviva	optiprint tray	optiprint zero	optiprint mock- up	optiprint match
<p>● Validierung abgeschlossen</p> <p>□ Validierung in Arbeit</p>												
Asiga Max (385 nm)	●	●	●	●	●	●	●	●	●	●	●	●
Asiga Pro 4K (385 nm)	●	●	●	●	●	●	●	●	●	●	●	●
EnvisionTEC D4K Pro (385 nm)	●	●	□	□	□	●	□	□	●	□	□	□
EnvisionTEC envisionOne (385 nm)	●	●	□	□	□	●	□	□	●	□	□	□
Ackuretta DENTIQ (405 nm)	●	●	●	●	●	●	●	●	●	●	●	●
Ackuretta SOL (405 nm)	●	●	●	●	●	●	●	●	●	●	●	●
Ackuretta FreeShape (405 nm)	●	●	●	●	●	●	●	●	●	●	●	●
HeyGears A2D (385 nm)	□	□	□	□	□	□	□	□	□	□	□	□
Microloy Eve (405 nm)	●	●	●	●	●	●	●	●	●	●	●	●
Microloy Versus (385 nm)	●	●	●	●	●	●	●	●	●	●	●	●
W2P SolLab 335/370/450 (385 nm)												
MiiCraft 125 (405 nm)	●		●	●					●			
MiiCraft Ultra (405 nm)												
Flashforge Hunter (405 nm)												
Phrozen Sonic XL 4K (405 nm)	●	●	●	●	□	●	□	□	●	□	□	□

Status: 01. Oktober 2021

Finden Sie regelmäßige updates unter www.dentona.de